

DESIGN

Veronique®
c o m p r e s s i o n w e a r

Compression Garments Overview

Version 08.03.1.



Why should I use a medical compression garment?

It is essential that the compression garment worn during the recovery period provides targeted compression, something that is not found in a standard girdle or athletic tights. Athletic tights, girdles and jog bras also do not have flat seams, an open crotch and are not intended for long-term wear.

Design Veronique® compression garments are very different from a pair of cycling shorts, a body girdle or bra you might see in your local store. Our garments are specifically designed for postoperative wear and have important features like flat, non-marking seams, inner panels that provide the targeted compression necessary during healing, and no labels or anything that will irritate your skin during extended wear.

About the company

Design Veronique® was founded over 22 years ago by Veronica C. Smith and is committed to designing, manufacturing and providing the finest quality compression garments to doctors, hospitals and patients throughout the world. We will not compromise on product quality or service. All garments are designed from the finest quality components and materials and manufactured in our California headquarters to maintain the highest manufacturing and quality control standards.

Product line

Design Veronique® has over 85 different styles of compression garments, available in three colors and a size range of XS through 4X for most styles. We also make custom garments if you don't find the style or size needed.

Finding the ideal garment for you

We make a lot of garments, we have the ideal garment to help you recover from any procedure. Which garment is best for you? Consult our garment finder at: www.designveronique.com/garmentfinder

The garment finder works by having you answer three questions: The type of procedure you had, the exact procedure performed and the number of days post-operative you are. Based on your responses, the garment finder will suggest the correct garment style and either a zipper or non-zippered style.



Ordering & Questions

Found the garment you need? You can order any of our products on our secure website: www.DesignVeronique.com or call 800.442.5800 and speak with one of our patient specialists. Please give us a call any time you have a question about one of our garments or a measurement.

Features of Zippered, First Stage Garments

Designed for immediate postoperative use

A zippered garment is worn during the first stage of recovery when you are swollen and tender. This phase typically ranges from immediately post-op to 10-14 days. The zippers allow you to change in and out of the garment with ease.

Sized to accommodate postoperative swelling and edema

The zippered style sizes are based on your measurements before surgery. Although you will be smaller after your procedure, postoperative swelling typically means you'll remain at your pre-op size for the first few weeks after your surgery. Using pre-op measurements ensures a proper, comfortable fit.

We recommend that you have more than one garment. Immediately post-op the doctor will put you in a zippered garment, this garment might become stained or soiled. In order to maintain continuous compression, the second garment can be worn while the first garment is being laundered.



#753
Zippered Above-Knee Girdle
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Lower abdominal, hip, and thigh procedures



#85
Zippered Torso Brief
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal procedures



#850
Zippered Abdominal Girdle
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal procedures



#B85
Zippered High-Back Torso Brief with Bra
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal, back & breast procedures



#B853
Zippered Above-Knee High-Back Girdle with Bra
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal, hip, thigh, back & breast procedures



#855
Zippered Full-Body Girdle
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal, hip, thigh, and knee procedures

Features of Non-Zippered, Second Stage Garments

Designed for body contouring and skin retracting

As your swelling begins to subside, a non-zippered garment will provide more targeted support and aid in skin retraction during the final body contouring phase.

Targeted Compression Panels

Non-zippered garments have strategically placed compression panels in key surgical areas to provide additional compression. The targeted compression provides more support and helps flush out any remaining fluids.

Designed for long-term wear


Non-zippered garments are particularly comfortable for long-term wear and allow you to resume regular daily activities.



#1453
Non-Zippered Above-Knee Girdle
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Lower abdominal, hip,
and thigh procedures



#165
Non-Zippered Torso Brief
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal procedures




#1650
Non-Zippered Abdominal Girdle
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal procedures



#1655
Non-Zippered Full-Body Girdle
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal, hip, thigh,
and knee procedures



#B165
Non-Zippered High-Back Torso
Brief with Bra
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal, back & breast
procedures



#B1653
Non-Zippered Above-Knee High-
Back Girdle with Bra
Sizes: XS - 4X
Colors: Black, White, Champagne
Ideal for: Abdominal, hip, thigh,
back & breast procedures

Breast Wear

Visit www.DesignVeronique.com to see all Breast Wear garments.




#450
D' Andrea — Optimum Support Cotton Knit Bra
Sizes: XS – 3X
Colors: Black, White
Ideal for: Breast reduction and reconstruction



#4510
Georgette—Fully Adjustable Implant Stabilizing Bra
Sizes: XS - 3X
Colors: Black, White, Champagne
Ideal for: Breast augmentation, stabilizing
and positioning implants



#454
Allyssandra — 2" Band Cotton Knit Bra
Sizes: XS – 3X
Colors: Black, White, Champagne
Ideal for: Breast augmentation and mastopexy



#510-X
Versatile Stabilizing Breast Wrap
Sizes: One Size
Colors: White
Ideal for: Breast augmentation, stabilizing,
and positioning implants



#457-Z
Yesmina-Z — Front Zippered Cotton Medical/Sports Bra
Sizes: XS – 3X
Colors: Black, White, Champagne
Ideal for: Breast augmentation, stabilizing and position-
ing implants



#510-J
3" Breast Wrap with Lace
Sizes: One Size
Colors: Black, White
Ideal for: Breast implant stabilization
and positioning

Upper Body Garments

Visit www.DesignVeronique.com to see all upper body garments.



#B943

Compression Arm Sleeve
with Adjustable Cotton Knit Bra

Sizes: XS - 4X

Colors: Black, White, Champagne

Ideal for: Arm and breast procedures



#BV933

Long Arm Sleeve Bolero

Sizes: XS - 4X

Colors: Black, White

Ideal for: Arm, back, and breast procedures



#350-9

9" Double Panel Binder

Sizes: XS - 4X

Colors: White

Ideal for: Abdominal procedures

Male Garments

Visit www.DesignVeronique.com to see all male garments.



#640

Adjustable Compression
Vest

Sizes: S - 4X

Colors: Black, White

Ideal for: Chest procedures



#1240

Compression Vest

Sizes: S - 4X

Colors: Black, White

Ideal for: Chest procedures

Facial Bands

Visit www.DesignVeronique.com to see all facial bands.



#210

Universal Facial Band[®]

Sizes: One Size

Colors: White, Champagne

Ideal for: Facial procedures



#210-N

Universal Facial and Neck Wrap

Sizes: M and L

Colors: White, Champagne

Ideal for: Facial and neck
procedures